

sempermed® 
classic

CLASSIC

THE ESTABLISHED CLASSIC

IDEAL TACTILE PROPERTIES – highly developed dipping process

HIGH SAFETY – ideal compounding

EXCELLENT GRIP – high comfort



THE ESTABLISHED CLASSIC

■ IDEAL TACTILE PROPERTIES

The highly developed dipping process guarantees constant wall thickness from fingers to cuff. This results in ideal tactile sensitivity and less hand fatigue even during longer operations.

■ SKIN CARE

The glove is powdered with pure cornstarch which makes it easy to don. The multiple washing process (leaching) makes the natural latex glove skin friendly and very comfortable to wear.

■ INCREASED SAFETY

Sempermed surgical gloves are classified as **Medical Device (MD)** as well as **Personal Protective Equipment (PPE)**. This means that the glove can also be used for handling chemicals¹⁾.

■ EXCELLENT GRIP AND COMFORT

Sensitivity and a safe grip are essential for a successful surgery. This glove has a micro-rough non-slip surface to provide a good grip when handling instruments.



TECHNICAL DATA

sempermed® classic

TYPE	Sterile surgical glove for single use, powdered		
MATERIAL	Natural rubber latex (NR)		
COLOUR	Natural white		
GLOVE SHAPE	Fully anatomical with rolled rim		
SIZE/OVERALL LENGTH as per EN 455-2	5 ½, 6 and 6 ½ 7, 7 ½ and 8 8 ½ and 9	min. 260 mm min. 270 mm min. 280 mm	
WALL THICKNESS measured in single layer in the palm area	min. 0.14 mm		
BARRIER PERFORMANCE Tightness as per EN 455-1	AQL 1.5		
FORCE AT BREAK as per EN 455-2	≥ 9 N		
DURABILITY	5 years in original packaging if stored as per DIN 7716, ISO 2230		
STERILISATION	Radiation STERILE R		
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223		
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Directive 89/686/EEC	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾		
PACKAGING	Left and right glove with turned up cuff in inner pouch, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 70 pairs In transport carton with sterilisation indicator: 420 pairs		
ARTICLE NUMBERS	Size 5 ½ 826054520 Size 6 826054600 Size 6 ½ 826054620 Size 7 826054700	Size 7 ½ 826054720 Size 8 826054800 Size 8 ½ 826054820 Size 9 826054900	



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
E-Mail: sempermed@semperitgroup.com
www.sempermed.com

Semperit Technische Produkte Gesellschaft m. b. H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
E-Mail: sempermed@semperitgroup.com · www.sempermed.com

¹⁾ Please note the Sempercare Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available under www.sempermed.com. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit does not accept any liability for incorrect use of the gloves. In case of doubt, obtain expert advice before use. The information and classification correspond to the latest status prior to printing. Technical details are average values from production and may vary in individual cases. Subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock.**